




PACT Antibody [1B9-1A7]

CATALOG NUMBER: 51-733



Immunohistochemistry of human tonsil tissue stained using PACT Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC, IP, WB
APPLICATIONS:	PACT antibody can be used in Western Blot starting at 1 - 3 ug/mL, and immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant protein
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	PACT antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	PRKRA, DYT16, HSD14, PKR-associating protein X, PACT, PKR-associated protein X
ACCESSION NO.:	O75569
PROTEIN GI NO.:	74735517
OFFICIAL SYMBOL:	PRKRA
GENE ID:	8575

Background

