



CA3 Antibody

CATALOG NUMBER: 70-267

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC
APPLICATIONS:	CA3 Antibody can be used for immunohistochemistry (tested on parafin embedded tissues) starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Recognizes Carbonic Anhydrase III (CAIII) .
IMMUNOGEN:	Recombinant Human CAIII.
HOST SPECIES:	Sheep

Properties

PURIFICATION:	Purified
PHYSICAL STATE:	Liquid
BUFFER:	20 mM PBS, 0.09% sodium azide, pH 7.2, 150 mM sodium chloride
STORAGE CONDITIONS:	CA3 antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Polyclonal
ISOTYPE:	IgG
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	CA3, Carbonate dehydratase III, Carbonic anhydrase 3, Carbonic anhydrase III, CA-III, CAIII, Car3
ACCESSION NO.:	P07451
PROTEIN GI NO.:	134047703
OFFICIAL SYMBOL:	CA3
GENE ID:	761

Background

FOR RESEARCH USE ONLY

December 13, 2016