



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

GPR151 Antibody

CATALOG NUMBER: 70-524

Specifications

SPECIES REACTIVITY:	Human, Mouse, Rat
TESTED APPLICATIONS:	ELISA, IF, IHC, WB
APPLICATIONS:	ELISA, IF, IHC-P (1:200), WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	GPR151 Antibody detects endogenous levels of total GPR151 protein.
IMMUNOGEN:	Synthetic peptide from human GPR151.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
BUFFER:	PBS (without Mg2+, Ca2+), pH 7.4, 150 mM sodium chloride, 0.02% sodium azide, 50% glycerol
STORAGE CONDITIONS:	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Polyclonal
ISOTYPE:	IgG
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	GPR151, AXOR40, G protein-coupled receptor 151, Galanin receptor 4, GALR4, GPCR-2037, HGPRBMY7, GALRL, PGR7
ACCESSION NO.:	Q8TDV0
PROTEIN GI NO.:	48428097
OFFICIAL SYMBOL:	GPR151
GENE ID:	134391

Background

FOR RESEARCH USE ONLY

December 13, 2016