



PELI1 Antibody

CATALOG NUMBER: 70-593

Specifications

SPECIES REACTIVITY:	Bat, Bovine, Chicken, Dog, Gibbon, Gorilla, Hamster, Horse, Human, Monkey, Mouse, Pig, Rabbit, Rat, Xenopus, Zebrafish
TESTED APPLICATIONS:	IHC, WB
APPLICATIONS:	PELI1 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 1:100, and Western Blotting starting at 1:500 - 1:1000.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	KLH-conjugated synthetic peptide corresponding to a portion of human Pellino 1 between amino acids 121-135.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Antiserum
BUFFER:	PBS containing 0.05% BSA and 0.05% sodium azide
STORAGE CONDITIONS:	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Polyclonal
ISOTYPE:	IgG
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	PELI1, Pellino 1, Pellino-1, Pellino (Drosophila) homolog 1, Pellino homolog 1 (Drosophila), PRISM, Protein pellino homolog 1, Pellino homolog 1
ACCESSION NO.:	Q96FA3
PROTEIN GI NO.:	37999756
OFFICIAL SYMBOL:	PELI1
GENE ID:	57162

Background

FOR RESEARCH USE ONLY

December 13, 2016