




FAP Antibody

CATALOG NUMBER: 48-343



Immunohistochemistry staining of FAP in uterus tissue using FAP Antibody.

Specifications

SPECIES REACTIVITY:	Gibbon, Gorilla, Human
TESTED APPLICATIONS:	IHC
APPLICATIONS:	FAP antibody can be used in Western Blot, immunohistochemistry starting at 10 ug/mL, and immunoprecipitation.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	FAP antibody was raised against a peptide located near the Extracellular domain of FAP (Human).
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.1% sodium azide.
STORAGE CONDITIONS:	FAP antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	FAP, DPPIV, FAP alpha, Fibroblast activation protein, FAPA, Seprase
ACCESSION NO.:	Q12884
PROTEIN GI NO.:	292495099
OFFICIAL SYMBOL:	FAP
GENE ID:	2191

Background

FOR RESEARCH USE ONLY

