




FOXO3 Antibody

CATALOG NUMBER: 48-953



Immunohistochemistry staining of FOXO3
in skin (formalin-fixed paraffin embedded)
tissue using FOXO3 Antibody.

Specifications

SPECIES REACTIVITY:	Bovine, Chicken, Dog, Human, Monkey, Mouse, Rat, Xenopus
TESTED APPLICATIONS:	IHC, WB
APPLICATIONS:	FOXO3 antibody can be used in ELISA, Western Blot, and immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	A portion of amino acids 600-673 of human FOXO3A.
IMMUNOGEN:	A portion of amino acids 600-673 of human FOXO3A.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.05% BSA, 0.05% sodium azide.
STORAGE CONDITIONS:	FOXO3 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	FOXO3, AF6q21, AF6q21 protein, Forkhead box O3, Forkhead box O3A, FOXO3A, FKHL1P2, Forkhead box protein O3, FOXO2, FKHL1
ACCESSION NO.:	O43524
PROTEIN GI NO.:	8134467
OFFICIAL SYMBOL:	FOXO3
GENE ID:	2309

Background

