





## PPP2R1B Antibody

CATALOG NUMBER: 49-619



Immunohistochemistry staining of PPP2R1B in cerebellum tissue using PPP2R1B Antibody.



Immunohistochemistry staining of PPP2R1B in kidney using PPP2R1B Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Bovine, Human, Mouse, Rat
<b>TESTED APPLICATIONS:</b>	IHC, WB
<b>APPLICATIONS:</b>	PPP2R1B antibody can be used in ELISA, Western Blot, and immunohistochemistry starting at 5 ug/mL.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	PPP2R1B antibody was raised against an internal peptide amino acids 7-19 of PPP2R1B (Rat).
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Ammonium Sulfate Precipitation
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, 0.08% sodium azide.
<b>STORAGE CONDITIONS:</b>	Store PPP2R1B antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	PPP2R1A, PP2A-Aalpha, PP2AAALPHA, PR65A
<b>ACCESSION NO.:</b>	P30153
<b>PROTEIN GI NO.:</b>	143811355
<b>OFFICIAL SYMBOL:</b>	PPP2R1A
<b>GENE ID:</b>	5518

### Background

FOR RESEARCH USE ONLY

