




ABCB7 Antibody

CATALOG NUMBER: 49-653



Immunohistochemistry staining of ABCB7
in thymus tissue using ABCB7 Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC
APPLICATIONS:	ABCB7 antibody can be used in Western Blot, and immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Amino acids 733 to 747 of human ABCB7
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Protein G Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.1% sodium azide.
STORAGE CONDITIONS:	ABCB7 antibody should be stored long term (months) at -20 °C and short term (weeks) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	ABCB7, ABC transporter 7 protein, ABC7, Atm1p, ATP-binding cassette 7, ASAT, EST140535
ACCESSION NO.:	O75027
PROTEIN GI NO.:	8928549
OFFICIAL SYMBOL:	ABCB7
GENE ID:	22

Background

FOR RESEARCH USE ONLY

