




PDGFRL Antibody

CATALOG NUMBER: 49-687



Immunohistochemistry staining of
PDGFRL in colon tissue using PDGFRL
Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC
APPLICATIONS:	PDGFRL antibody can be used in Western Blot, and immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Amino acids 343 to 357 of human PDGFRL
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Protein G Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.09% sodium azide.
STORAGE CONDITIONS:	PDGFRL antibody should be stored long term (months) at -20 °C and short term (weeks) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	PDGFRL, PDGRL, PRLTS, PDGFR-like protein
ACCESSION NO.:	Q15198
PROTEIN GI NO.:	74762141
OFFICIAL SYMBOL:	PDGFRL
GENE ID:	5157

Background

FOR RESEARCH USE ONLY

