




## T-box 4 Antibody

CATALOG NUMBER: 49-740



Immunohistochemistry staining of T-box 4 in skeletal muscle (formalin-fixed paraffin embedded) tissue using T-box 4 Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	ELISA, IHC
<b>APPLICATIONS:</b>	T-box 4 antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 ug/mL, and immunoprecipitation
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	Amino acids 529 to 545 of human TBX4
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Protein G Column
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, 0.09% sodium azide.
<b>STORAGE CONDITIONS:</b>	T-box 4 antibody should be stored long term (months) at -20 °C and short term (weeks) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	TBX4, SPS, T-box 4, T-box protein 4
<b>ACCESSION NO.:</b>	P57082
<b>PROTEIN GI NO.:</b>	51338786
<b>OFFICIAL SYMBOL:</b>	TBX4
<b>GENE ID:</b>	9496

### Background

FOR RESEARCH USE ONLY

