



LYPD3 Antibody

CATALOG NUMBER: 49-988

Specifications

SPECIES REACTIVITY:	Gibbon, Gorilla, Human
TESTED APPLICATIONS:	IHC
APPLICATIONS:	LYPD3 antibody can be used in Western Blot, immunohistochemistry starting at 10 ug/mL, immunoprecipitation, and flow cytometry starting at 1:25 - 1:50.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	BLAST analysis of the peptide immunogen showed no homology with other human proteins.
IMMUNOGEN:	LYPD3 antibody was raised against C-Terminal domain of human (Human).
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.1% sodium azide.
STORAGE CONDITIONS:	LYPD3 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	LYPD3, C4.4A, MIG-C4, LY6/PLAUR domain containing 3
ACCESSION NO.:	O95274
PROTEIN GI NO.:	90185271
OFFICIAL SYMBOL:	LYPD3
GENE ID:	27076

Background

FOR RESEARCH USE ONLY

December 13, 2016