




YY1 Antibody [4D2]

CATALOG NUMBER: 51-820



Immunohistochemistry of human testis tissue stained using YY1 Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IF, IHC, WB
APPLICATIONS:	YY1 antibody can be used in Western Blot starting at 1 - 3 ug/mL, immunohistochemistry starting at 2.5 ug/mL, and immunoprecipitation.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant protein
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	YY1 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	YY1, Delta transcription factor, INO80S, NF-E1, INO80 complex subunit S, YY1 transcription factor, YIN-YANG-1, YY-1, Yin and yang 1, UCRBP, Yin and Yang 1 protein
ACCESSION NO.:	P25490
PROTEIN GI NO.:	3915889
OFFICIAL SYMBOL:	YY1
GENE ID:	7528

FOR RESEARCH USE ONLY

December 13, 2016