



## ABCD1 Antibody

CATALOG NUMBER: 70-396

### Specifications

<b>SPECIES REACTIVITY:</b>	Human, Mouse
<b>TESTED APPLICATIONS:</b>	IHC, WB
<b>APPLICATIONS:</b>	ABCD1 Antibody can be used for Immunohistochemistry in paraffin embeded tissues starting at 5 ug/mL, and Western Blotting.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	Recombinant protein
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Affinity Purified
<b>BUFFER:</b>	PBS, pH 7.4, 0.02% sodium azide
<b>CONCENTRATION:</b>	1 mg/ml
<b>STORAGE CONDITIONS:</b>	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
<b>CLONALITY:</b>	Polyclonal
<b>ISOTYPE:</b>	IgG
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	ABCD1, ABC42, Adrenoleukodystrophy protein, ALD, Aldgh, Adrenoleukodystrophy, ALDP, X-ALD
<b>ACCESSION NO.:</b>	P33897
<b>PROTEIN GI NO.:</b>	67476960
<b>OFFICIAL SYMBOL:</b>	ABCD1
<b>GENE ID:</b>	215

### Background

FOR RESEARCH USE ONLY

December 13, 2016