



## LAMB1 Antibody

CATALOG NUMBER: 70-405

### Specifications

#### SPECIES REACTIVITY:

**TESTED APPLICATIONS:** IHC, WB

**APPLICATIONS:** ELISA, IF, IHC-P (1:100), WB

**USER NOTE:** Optimal dilutions for each application to be determined by the researcher.

**SPECIFICITY:** LAMB1 Antibody detects endogenous levels of total LAMB1 protein.

**IMMUNOGEN:** Synthetic peptide from human LAMB1.

**HOST SPECIES:** Rabbit

### Properties

**PURIFICATION:** Immunoaffinity Chromatography

**BUFFER:** PBS (without Mg2+, Ca2+), pH 7.4, 150 mM sodium chloride, 0.02% sodium azide, 50% glycerol

**STORAGE CONDITIONS:** Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.

**CLONALITY:** Polyclonal

**ISOTYPE:** IgG

**CONJUGATE:** Unconjugated

### Additional Info

**ALTERNATE NAMES:** LAMB1, CLM, Laminin B1 chain, Laminin, beta 1, Laminin-12 subunit beta, Laminin-6 subunit beta, Laminin-1 subunit beta, Laminin-8 subunit beta, Laminin beta-1, Laminin beta-1 chain, Laminin subunit beta-1, Laminin-10 subunit beta, Laminin-2 subunit beta

**ACCESSION NO.:** P07942

**PROTEIN GI NO.:** 317373377

**OFFICIAL SYMBOL:** LAMB1

**GENE ID:** 3912

### Background

FOR RESEARCH USE ONLY

December 13, 2016