



FOXD3 Antibody

CATALOG NUMBER: 70-406

Specifications

| | |
|-----------------------------|--|
| SPECIES REACTIVITY: | Human, Mouse |
| TESTED APPLICATIONS: | ELISA, IHC, WB |
| APPLICATIONS: | ELISA, IHC-P (1:100), WB |
| USER NOTE: | Optimal dilutions for each application to be determined by the researcher. |
| SPECIFICITY: | FOXD3 Antibody detects endogenous levels of total FOXD3 protein. |
| IMMUNOGEN: | Synthetic peptide from human FOXD3. |
| HOST SPECIES: | Rabbit |

Properties

| | |
|----------------------------|---|
| PURIFICATION: | Immunoaffinity Chromatography |
| BUFFER: | PBS (without Mg2+, Ca2+), pH 7.4, 150 mM sodium chloride, 0.02% sodium azide, 50% glycerol |
| STORAGE CONDITIONS: | Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles. |
| CLONALITY: | Polyclonal |
| ISOTYPE: | IgG |
| CONJUGATE: | Unconjugated |

Additional Info

| | |
|-------------------------|---|
| ALTERNATE NAMES: | FOXD3, AIS1, Drosophila forkhead homologue, Genesis, Forkhead box protein D3, Forkhead box D3, HFH2, VAMAS2 |
| ACCESSION NO.: | Q9UJU5 |
| PROTEIN GI NO.: | 8134475 |
| OFFICIAL SYMBOL: | FOXD3 |
| GENE ID: | 27022 |

Background

FOR RESEARCH USE ONLY

December 13, 2016