



CSF2RB / CD131 Antibody

CATALOG NUMBER: 70-451

Specifications

SPECIES REACTIVITY:	Human, Mouse, Rat
TESTED APPLICATIONS:	ELISA, IHC, WB
APPLICATIONS:	ELISA, IHC-P (1:100), WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	IL-3R beta antibody detects endogenous levels of total IL-3R beta protein.
IMMUNOGEN:	Synthetic peptide from human IL-3R beta.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
BUFFER:	PBS (without Mg2+, Ca2+), pH 7.4, 150 mM sodium chloride, 0.02% sodium azide, 50% glycerol
STORAGE CONDITIONS:	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Polyclonal
ISOTYPE:	IgG
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	CSF2RB, CD131 antigen, CDw131, IL3RB, IL5RB, SMDP5, CD131, GM-CSF receptor beta chain
ACCESSION NO.:	P32927
PROTEIN GI NO.:	1345923
OFFICIAL SYMBOL:	CSF2RB
GENE ID:	1439

Background

FOR RESEARCH USE ONLY

December 13, 2016