



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

K-562 Cell Slide

CATALOG NUMBER: 17-004

Specifications

TISSUE CELL TYPE:	K-562
SPECIES:	Human
GENERAL DIAGNOSIS:	K-562 Cell Slide
AVAILABLE LOTS:	human (bone marrow) chronic myelogenous leukemia
TESTED APPLICATIONS:	ICC
APPLICATIONS:	K-562 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

Properties

FIXATION:	4% paraformaldehyde
STORAGE CONDITIONS:	K-562 cell slide can be stored for up to 6 months at 4°C under desiccate conditions.

Background

K-562 (bone marrow; chronic myelogenous leukemia)

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

FOR RESEARCH USE ONLY

December 12, 2016