



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

RAW 264.7 Cell Slide

CATALOG NUMBER: 17-203

Specifications

TISSUE CELL TYPE:	RAW 264.7
SPECIES:	Mouse
GENERAL DIAGNOSIS:	RAW 264.7 Cell Slide
AVAILABLE LOTS:	BALB/c Mouse (adult, Male) ascites, Abelson murine leukemia virus transformed, macrophage
TESTED APPLICATIONS:	ICC
APPLICATIONS:	RAW 264.7 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

Properties

FIXATION:	4% paraformaldehyde
STORAGE CONDITIONS:	RAW 264.7 cell slide can be stored for up to 6 months at 4°C under desiccated conditions.

Background

RAW 264.7 (Abelson murine leukemia virus-induced tumor)

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

FOR RESEARCH USE ONLY

December 12, 2016