



CD83 Recombinant Protein

CATALOG NUMBER: 92-645

Specifications

SPECIES:	Mouse
SOURCE SPECIES:	Human Cells
SEQUENCE:	Met22-Ala134
FUSION TAG:	C-6 His tag
APPLICATIONS:	This recombinant protein can be used for biological assays. For research use only.

Properties

PURITY:	Greater than 95% as determined by reducing SDS-PAGE. Endotoxin level less than 0.1 ng/ug (1 IEU/ug) as determined by LAL test.
PREDICTED MOLECULAR WEIGHT:	17.4 kD
PHYSICAL STATE:	Lyophilized
BUFFER:	Lyophilized from a 0.2 um filtered solution of PBS, pH7.4. It is not recommended to reconstitute to a concentration less than 100 ug/ml. Dissolve the lyophilized protein in ddH ₂ O.
STORAGE CONDITIONS:	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Additional Info

ALTERNATE NAMES:	CD83 Antigen, hCD83, B-Cell Activation Protein, Cell Surface Protein HB15, CD83, HB15
ACCESSION NO.:	O88324

Background

CD83 antigen is a single-pass type I membrane protein which contains one Ig-like V-type (immunoglobulin-like) domain. CD83 is expressed by activated lymphocytes, Langerhans cells and interdigitating reticulum cells. It contains one Ig-like V-type (immunoglobulin-like) domain. , the soluble CD83 has the opposite effect and has an immune inhibitory capacity. Due to its immune inhibitory function, CD83 may play a significant role in antigen presentation or the cellular interactions that follow lymphocyte activation.

FOR RESEARCH USE ONLY

December 14, 2016