



Interior

Front View



Side View



**Cat# 10819-780/HCLP-10S**

115VAC/60Hz/1PH  
NEMA-5-15P

| Cu.Ft. | Defrost | Ft Door | Int Door | Shelves                 | Interior Dimensions |        |        | Refrigerant | H.P. | Amps | Weight   |
|--------|---------|---------|----------|-------------------------|---------------------|--------|--------|-------------|------|------|----------|
|        |         |         |          |                         | W"                  | D"     | H"     |             |      |      |          |
| 10.5   | Cycle   | 1 Solid | 0        | 6 Adj<br>+Footed Bottom | 19 7/8              | 18 3/8 | 53 1/2 | R600a       | 1/3  | 2.3  | 145 lbs. |

Note: Interior dimensions are estimates and subject to change without notice