



Interior

Front View



Side View



Cat# 10819-540/HCLP-69-TS

115VAC/60Hz/1PH
NEMA-5-15P

Cu.Ft.	Defrost	Door	Int Door	Shelves	Interior Dimensions			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
69	Cycle	3 Slide	0	12 Adj	74 1/4	21 3/4	56 1/2	R290	1/2	9.3	645 lbs.

Note: Interior dimensions are estimates and subject to change without notice