



Interior

Front View



Side View



Cat# 10820-380/SLS-12

115VAC/60Hz/1PH
NEMA-5-15P

Cu.Ft.	Defrost	Door	Int Door	Shelves	Interior Dimensions			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
12	Cycle	1 Solid	0	3 Adj	19 ³ / ₄	19 ¹ / ₄	36 ¹ / ₄	R134a	³ / ₈	6	187 lbs.

Note: Interior dimensions are estimates and subject to change without notice