



307PR

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name

OneQuant™ BSA

(Immunological Grade)

(Cat. # 786-090)



think proteins! think G-Biosciences www.GBiosciences.com

INTRODUCTION

Immunological grade BSA is a high purity grade bovine serum albumin (BSA) prepared for immunotechnology research applications such as a carrier protein for antibody production. Immunogenic peptides are coupled to BSA with the use of protein cross-linkers, normally to primary amine residues on BSA. BSA (MW: 67kD) has a total of 59 lysine amine groups of which ~30-35 are sterically available for coupling. BSA is supplied lyophilized in an *OneQuant*[™] format.

ITEM(S) INCLUDED

Description	Size
OneQuant [™] BSA	8 x 2mg/vial

STORAGE CONDITIONS


It is shipped at ambient temperature. Upon receiving store at 4°C.

PROTOCOL

Reconstitute with pure water or a buffer of choice to desired concentration.

RELATED PRODUCTS

Download our Antibody Production Handbook.



<http://info.gbiosciences.com/complete-Antibody-Production-handbook>

For other related products, visit our website at www.GBiosciences.com or contact us.

Last saved: 2/14/2013 CMH

This page is intentionally left blank



www.GBiosciences.com