



AffinityPrep Resins

ImmunoReagents' Antibody Cross-Adsorption Resins are used to reduce assay interference caused from anti-animal (AA) and heterophilic antibodies. These resins are routinely used when multiple antibodies (i.e., sandwich ELISA) are used in an assay, and antibody-antibody interactions are undesirable. Each resin is prepared using IR's exclusive AffinityPrep technology incorporating highly purified species specific immunoglobulins. Kits provide all the necessary reagents and materials for fast, convenient cross absorption of your samples. Kits include: buffers, reagents, chromatography column, centrifuge tubes, and instructions.



Antibody Cross-Absorption Resins and Kits

Product	Resin Cat. No.	Kit Cat. No.
Bovine IgG-AffinityPrep	10148-968	10018-804
Chicken IgY-AffinityPrep	10148-970	10018-690
Donkey IgG-AffinityPrep	10148-972	10018-692
Goat IgG-AffinityPrep	10148-974	10018-694
Horse Serum-AffinityPrep	10013-954	10018-696
Human IgG-AffinityPrep	10148-976	10018-698
Llama IgG-AffinityPrep	10148-978	10018-702
Mouse IgG-AffinityPrep	10148-980	10018-704
Rabbit IgG-AffinityPrep	10148-982	75806-250
Sheep IgG-AffinityPrep	10148-984	10018-708

PRIMARY ANTIBODIES

SECONDARY ANTIBODIES

ANTIBODY PURIFICATION REAGENTS



ImmunoReagents, Inc.

BULK IMMUNOASSAY REAGENTS

ELISA & WESTERN BLOT REAGENTS

PURIFIED PROTEINS & IMMUNOGLOBULINS



1.800.932.5000
vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

0617 Lit. No. 090001W