

Think. Inside. The. Box.

PURELAB® Ultra Bioscience*



DI pack

LC182/LC183

Lifetime expectancy: 6 months

LC184/LC185

What is purified from the water?

Ions

Impact of not changing:

Weakly bound ions are released back into the water reducing resistivity



In-line Filtration

LC151

Lifetime expectancy: 24 months

What is purified from the water?

Bacteria and microbial by-products such as endotoxins

Impact of not changing:

Bacterial and endotoxin specification will not be maintained

Reduced flow rate



Point of Use Filter**

LC134

Lifetime expectancy: 6 months

What is purified from the water?

Particulates and bacteria

Impact of not changing:

Potential recontamination of purified water and reduced flow rates

*A combination of treatments and technologies produce the final water specification

** For DNase, RNase & endotoxin specification at point of use, please use a LC197 Biofilter