



Product Data Sheet

Product Name: Substance P, FAM-labeled

Catalog Number: AS-61201 (1 mg) **Lot Number:** See label on vial

Sequence: FAM-Arg-Pro-Lys-Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-NH₂ (3-letter code)
FAM-RPKPQQFFGLM-NH₂ (1-letter code)

Molecular Weight: 1706.0

Peptide Purity: >95%

Appearance: Lyophilized yellow powder

Peptide Reconstitution: Substance P, FAM-labeled peptide is freely soluble in water.

Storage: Substance P, FAM-labeled peptide is shipped at ambient temperature. Upon receipt, store lyophilized peptide at –20°C or lower. Reconstituted peptide can be aliquoted and stored at –20 °C or lower.

Description: A fluorescent (FAM)-labeled Substance P, Abs/Em = 494/521 nm.

Related Products:

Name	Cat #	Size
Substance P	AS-24279	1 mg
RPKPQQFFGLM-NH ₂	AS-24280	5 mg

Substance P (5-11), amide QQFFGLM-NH ₂	61692	1 mg
Substance P (7-11) FFGLM-NH ₂	61693	1 mg
Substance P (8-11) FGLM-NH ₂	61694	1 mg
Substance P (9-11) GLM-NH ₂	61695	1 mg

For Research Use Only