

The New Boekel Multi-Purpose Tube Rotator

- Digital microprocesor control
- Adjustable speed provides gentle or vigorous mixing action
- Adjustable angled head offers mixing angles from 90 to 180 degrees
- Five carousel options provide greater versatility for the most tube mixing needs
- Timer features offers time-down or time-up counting feature
- 2 year warranty
- Power pac is UL, CSA and CE approved
- Includes 260653 - 14 Place Tube Carousel for 50mL



Catalog No.:	260650
Speed Range:	4 - 80 rpm
Carousel angle	90 - 180 degrees
Timer (with auto shut off)	Up to 99 hrs and 59 minutes
Carousel capacity	(46) 9 - 11 mm tubes, (28) 14 - 18 mm tubes, (14) 27 - 34 mm tubes, (63) 18.5 mm culture tubes (drum rotor) and (63) 14 mm culture tubes (drum rotor)
Electrical:	19.4V dc, 474A via 100-240 AC adapter
Accessories	
260651	46 place tube carousel for 1.5/ 2.0 ml tubes
260652	28 place tube carousel for 15 mL tubes
260653	14 place tube carousel for 50 mL tubes
260654	63 place tube carousel for 18.5 mm cul. tubes
260655	63 place tube carousel for 14 mm cul. tubes