

**MATERIAL SAFETY DATA SHEET
TRYPTIC SOY AGAR MEDIA**

Page 1 of 3

1. Product And Company Identification

Teknova Inc.

2290 Bert Dr.

Hollister, Ca 95023

Telephone Number: 831-637-1100

Emergency Contacts/Telephone Numbers

Chemtrec Emergency Phone Number: (800) 424-9300

Issue Date: 7/30/2009

Catalog Number: T0520

Product Name: TRYPTIC SOY AGAR MEDIA

2. Composition/Information On Ingredients

Ingredients	Cas Number	Sara #:
TRYPTIC SOY AGAR	N/A	N/A

3. Hazards Identification

Hmis Rating

Health: 1

Flammability: 0

Reactivity: 0

Nfpa Rating:

Health: 1

Flammability: 0

Reactivity: 0

Refer To Section 11 For More Information.

4. First Aid Measures

Eye

In Case Of Contact, Immediately Flush Eyes With Water. Assure Adequate Flushing By Separating The Eyelids With Fingers.

Skin

In Case Of Contact, Wash Skin With Soap And Water. Wash Contaminated Clothing Before Reuse.

Inhalation

If Inhaled, Remove To Fresh Air. If Breathing Is Difficult, Give Oxygen.

Ingestion

If Swallowed, Wash Out Mouth With Water Provided Person Is Conscious.

5. Firefighting Measures

Extinguishing Media

Water Spray. Carbon Dioxide, Dry Chemical Powder Or Appropriate Foam.

Firefighting Procedures

Wear Self-Contained Breathing Apparatus And Protective Clothing To Prevent Contact With Skin And Eyes.

6. Accidental Release Measures

Personal Precaution

Avoid forming dust.

Environmental precautions

Do not let product enter drains.

Methods For Cleaning Up

Sweep Up, Keep In A Closed Container And Hold For Waste Disposal.

7. Handling And Storage

Handling

Ventilate areas where dust is formed. Follow normal measure to prevent fire.

Storage:

Suitable: Keep Tightly Closed in a dry area that is well ventilated.

**MATERIAL SAFETY DATA SHEET
TRYPTIC SOY AGAR MEDIA**

Page 2 of 3

Recommended storage temperature: 2 - 8 °C
hygroscopic

8. Exposure Controls/Personal Protection

Personal Protective Equipment:

Wear Appropriate Niosh/Msha-Approved Respirator, Chemical-Resistant Gloves, Safety Goggles, Other Protective Clothing.

General Hygiene Measures:

Wash Thoroughly After Handling.

9. Physical And Chemical Properties

Appearance And Physical Properties: Solid

N/A

10. Stability And Reactivity

Stability

Stable: Stable.

Hazardous Decomposition Products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Hydrogen chloride gas, Sodium Oxides.

11. Toxicological Information

Signs And Symptoms Of Exposure

To The Best Of Our Knowledge, The Chemical, Physical, And Toxicological Properties Have Not Been Thoroughly Investigated.

Potential Health Effects

Inhalation May be harmful if inhaled. May cause respiratory tract irritation.

Skin May be harmful if absorbed through skin. May cause skin irritation.

Eyes May cause eye irritation.

Ingestion May be harmful if swallowed.

12. Ecological Information

Data Not Yet Available.

13. Disposal Considerations

Dispose of any unused product. Observe All Federal, State And Local Environmental Regulations.

14. Transport Information

Dot

Proper Shipping Name: None

Non-Hazardous For Transport: This Substance Is Considered To Be Non-Hazardous For Transport.

Iata

Non-Hazardous For Air Transport: Non-Hazardous For Air.

15. Regulatory Information

United States Regulatory Information

Sara Listed: No

Canada Regulatory Information

Whmis Classification: This Product Has Been Classified In Accordance With The Hazard Criteria Of The Cpr, And The Msds Contains All The Information Required By The Cpr.

Dsl: No

Ndsl: No

16. Other Information

Disclaimer

For R&D Use Only. Not For Drug, Household Or Other Uses.

Warranty

The Above Information Is Believed To Be Correct But Does Not Purport To Be All Inclusive And Shall Be Used Only As A Guide. Teknova, Inc. Shall Not Be Held Liable For Any Damage Resulting From Handling Or From Contact With The Above Product.

MATERIAL SAFETY DATA SHEET
TRYPTIC SOY AGAR MEDIA

Page 3 of 3

Teknova, Inc.
