

Name of the Product : **Lactobacillus Selection Oxgall Agar Base**  
 Code No. : **M1165**  
 Section 1 : **Chemical Identification**  
     Code No. : M1165  
     Name of the Product : Lactobacillus Selection Oxgall Agar Base  
     Produced by : HiMedia Laboratories Pvt. Ltd.  
     Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India.  
     Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468

## Section 2 : Composition

Ingredients	Grams/Litre	Ingredients	Grams/Litre
Casein enzymic hydrolysate	10.00	Yeast extract	5.00
Sodium acetate	25.00	Dextrose	20.00
Monopotassium hydrogen phosphate	6.00	*Ammonium citrate	2.00
Oxgall	1.50	Polysorbate 80	1.00
Magnesium sulphate	0.575	*Ferrous sulphate	0.034
*Manganese sulphate	0.12	Agar	15.00

Final pH (at 25°C) 5.5±0.2

\* Irritant / Harmful material



Ammonium citrate

CAS No. : 3458-72-8

R : 36/37/38 S : 26-36



Ferrous sulphate

CAS No. : 7782-63-0

R : 22-36/37/38 S : 26-36

RTECS : NO 8510000



Manganese sulphate

CAS No. : 10034-96-5

R : 48/20/22 S : 22

RTECS : OP 0893500

## Section 3 : Hazards Identification



Hazard : This symbol designates substances which may have an irritant effect on skin, eyes and respiratory organs.

Caution : Do not breathe vapours and avoid contact with skin and eyes.



Hazard : Inhalation and ingestion of, or skin penetration by these substances is harmful to one's health. Non-recurring, recurring or lengthy exposure to these substances may result in irreversible damage.

Caution : Avoid contact with the human body, including inhalation of the vapours and in cases of malaise consult a doctor.

---

**Section 4 : First - Aid Measures**

In case of swallowing, wash out mouth with water provided person is conscious. In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. In case of inhalation, remove to fresh air. If breathing becomes difficult, call physician. In case of contact with eyes, wash with copious amounts of water for at least 15 minutes. Assure adequate washing by separating the eyelids with fingers. In serious conditions, call a physician, show the container or label.

---

**Section 5 : Fire Fighting Measures**

Extinguishing Media. Water spray.

Carbon dioxide, dry chemical powder or appropriate foam. Special fire-fighting procedures. Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

---

**Section 6 : Accidental Release Measures**

Wear respirator, chemical safety goggles, rubber boots and heavy rubber gloves. Sweep up, place in a bag and hold for waste disposal. Avoid raising dust. Ventilate area and wash spill site after material pickup is complete.

---

**Section 7 : Handling and Storage**

Handling - Refer to Section 8

Storage - Store below 8°C

---

**Section 8 : Exposure Controls / Personal Protection**

Wear appropriate Niosh/Msha-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.

---

**Section 9 : Physical and Chemical Properties**

Appearance : Homogeneous powder

Colour : Yellow

---

**Section 10 : Stability and Reactivity**

Stability : Product is stable if stored as per the conditions specified under storage of Section No. 7. Product loses its potency/performance above 45°C.

Conditions to avoid : Heat and light.

Hazardous polymerization will not occur.

---

**Section 11 : Toxicological Information**

Acute Effects : May be harmful if swallowed.

Exposure can cause : Stomach pains, vomiting, diarrhoea. Prolonged or repeated exposure may cause allergic reactions in certain sensitive individuals.

RTECS No. : NO 8510000 OP 0893500

For complete information RTECS may be referred.

---

**Section 12 : Ecological Information**

Data not available

---

Section 13 : **Disposal Considerations**

On completion of work all used or unusable preparations of this product and derivatives thereof are to be disposed off by autoclaving and/or by incineration. Dispose off waste in accordance with all applicable Federal, State and local laws.

---

Section 14 : **Transport Information**

UN No. : Not applicable.

---

Section 15 : **Regulatory Information**

Risk Phrases : Irritating to eyes, respiratory system and skin

Risk Phrases : Harmful if swallowed

Risk Phrases : Harmful : danger of serious damage to health by prolonged exposure through inhalation and if swallowed

Safety Phrases: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

Safety Phrases: Wear suitable protective clothing

Safety Phrases: Do not breathe dust



Significance of signs : This symbol designates substances which may have an irritant effect on skin, eyes and respiratory organs.



Significance of signs : Inhalation and ingestion of, or skin penetration by these substances is harmful to one's health. Non-recurring, recurring or lengthy exposure to these substances may result in irreversible damage.

---

Section 16 : **Other Information**

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for user's obligation to investigate and determine the suitability of the information for their particular purpose.