

Name of the Product : **FAGI Agar**  
 Code No. : **M1196**  
 Section 1 : **Chemical Identification**  
     Code No. : M1196  
     Name of the Product : FAGI Agar  
     Produced by : HiMedia Laboratories Pvt. Ltd.  
     Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India.  
     Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468

Section 2 : **Composition**

Ingredients	Grams/Litre	Ingredients	Grams/Litre
Peptic digest of animal tissue	3.00	Magnesium sulphate	0.20
*Manganese sulphate	0.05	Agar	20.00
Final pH (at 25°C) 6.9±0.2			

\* Harmful / Dangerous for Environment material



Manganese sulphate  
 CAS No. : 10034-96-5  
 R : 48/20/22-51/53 S : 22-61  
 RTECS : OP 0893500

Section 3 : **Hazards Identification**

Hazard : Inhalation and ingestion of, or skin penetration by these substances is harmful to one's health. Non-recurring, recurring or lengthy exposure to these substances may result in irreversible damage.  
 Caution : Avoid contact with the human body, including inhalation of the vapours and in cases of malaise consult a doctor.

Hazard : Substances which are harmful to the aquatic (aquatic organisms, waters), as well as the non-aquatic environment (fauna, flora, atmosphere) or which have a detrimental effect at longer term.  
 Caution : Avoid release to the environment.

Section 4 : **First - Aid Measures**

In case of swallowing, wash out mouth with water provided person is conscious. In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. In case of inhalation, remove to fresh air. If breathing becomes difficult, call physician. In case of contact with eyes, wash with copious amounts of water for at least 15 minutes. Assure adequate washing by separating the eyelids with fingers. In serious conditions, call a physician, show the container or label.

Section 5 : **Fire Fighting Measures**

Extinguishing Media. Water spray.  
 Carbon dioxide, dry chemical powder or appropriate foam. Special fire-fighting procedures. Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

Section 6 : **Accidental Release Measures**

Wear respirator, chemical safety goggles, rubber boots and heavy rubber gloves. Sweep up, place in a bag and hold for waste disposal. Avoid raising dust. Ventilate area and wash spill site after material pickup is complete.

Section 7	<b>: Handling and Storage</b> Handling - Refer to Section 8 Storage - Store below 30°C
Section 8	<b>: Exposure Controls / Personal Protection</b> Wear appropriate NIOSH/MSHA-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.
Section 9	<b>: Physical and Chemical Properties</b> Appearance : Homogeneous powder Colour : Light Yellow
Section 10	<b>: Stability and Reactivity</b> Stability : Product is stable if stored as per the conditions specified under storage of Section No. 7. Product loses its potency/performance above 45°C. Conditions to avoid : Heat and light. Hazardous polymerization will not occur.
Section 11	<b>: Toxicological Information</b> Acute Effects : May be harmful if swallowed. Exposure can cause : Stomach pains, vomiting, diarrhoea. Prolonged or repeated exposure may cause allergic reactions in certain sensitive individuals. RTECS No. : OP 0893500 For complete information RTECS may be referred.
Section 12	<b>: Ecological Information</b> Data not available
Section 13	<b>: Disposal Considerations</b> On completion of work all used or unusable preparations of this product and derivatives thereof are to be disposed off by autoclaving and/or by incineration. Dispose off waste in accordance with all applicable Federal, State and local laws.
Section 14	<b>: Transport Information</b> UN No. : Not applicable.
Section 15	<b>: Regulatory Information</b> Risk Phrases : Harmful : danger of serious damage to health by prolonged exposure through inhalation and if swallowed Risk Phrases : Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment Safety Phrases : Do not breathe dust Safety Phrases: Avoid release to the environment. Refer to special instructions/Safety data sheets Significance of signs : Inhalation and ingestion of, or skin penetration by these substances is harmful to one's health. Non-recurring, recurring or lengthy exposure to these substances may result in irreversible damage. Significance of signs : Substances which are harmful to the aquatic (aquatic organisms, waters), as well as the non-aquatic environment (fauna, flora, atmosphere) or which have a detrimental effect at longer term.





Section 16 : **Other Information**

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for user's obligation to investigate and determine the suitability of the information for their particular purpose.

---

•