



Item Number BR132
Item Bromocresol Green, Reagent, ACS
CAS Number 76-60-8
Molecular Formula $C_{21}H_{14}Br_4O_5S$
Molecular Weight 698.02
Clarity of Solution..... To pass test
Visual Transition Interval..... pH 3.8 (yellow) to pH 5.4 (blue)

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: 14422 S. San Pedro St.
Gardena, CA 90248
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843

East Coast Facility: 755 Jersey Ave.
New Brunswick, NJ 08901