



Item Number C1155
Item Calcium Phosphate Tribasic, Powder, Reagent
CAS Number 7758-87-4
Molecular Formula $\text{Ca}_{10}(\text{OH})_2(\text{PO}_4)_6$
Molecular Weight 1004.6

MAXIMUM LIMITS

Arsenic (As).....	0.0001%
Chloride (Cl).....	0.05%
Lead (Pb).....	2 mg/kg
Insoluble in Hydrochloric Acid.....	0.07%
Iron (Fe).....	0.03%
Nitrogen (as N).....	0.05%
Sulfate (SO_4).....	0.3%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843