



Item Number C1228
Item Choline Bitartrate, FCC
CAS Number 87-67-2
Molecular Formula C₉H₁₉NO₇
Molecular Weight 253.25
Assay (C₉H₁₉NO₇; Anhydrous Basis)..... Min. 98.0%
Optical Rotation..... +17.5 to +18.5°
MAXIMUM LIMITS
1,4-Dioxane..... To pass test
Lead..... 2 mg/kg
Residue on Ignition..... 0.1%
Water..... 0.5%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843