



Item Number	C1388
Item	Calcium Peroxide, FCC
CAS Number	1305-79-9
Molecular Formula	CaO ₂
Molecular Weight	72.08
Assay (CaO ₂)	60.0%
MAXIMUM LIMITS	
Fluoride (F)	0.005%
Lead (Pb)	4 mg/kg

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843