



Item Number C1403
Item Creatinine, NF
CAS Number 60-27-5
Molecular Formula $C_4H_7N_3O$
Molecular Weight 113.12
Assay (Dried basis)..... 98.5-102.0%
MAXIMUM LIMITS
Loss on Drying..... 3.0%
Residue on Ignition..... 0.2%
Heavy Metals..... 0.001%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843