



Item Number C1479

Item Carbomer 1342, NF

CAS Number 9003-01-4

Molecular Formula

Molecular Weight

Carboxylic Acid Groups (Dried Basis).....52.0-62.0%

Viscosity of a Neutralized 1.0% Aqueous Dispersion..... 9,500-26,500 cP

MAXIMUM LIMITS

Benzene..... 0.2%

Loss on Drying..... 2.0%

Heavy Metals.....0.002%

Organic Volatile Impurities..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843