



Item Number C1659
Item Chlorambucil, USP
CAS Number 305-03-3
Molecular Formula $C_{14}H_{19}Cl_2NO_2$
Molecular Weight 304.21

Assay (Anhydrous Basis).....98.0-101.0%
Melting Range..... 65°-69°C

MAXIMUM LIMITS

Water..... 0.5%
Residual Solvents.....To pass test

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	