



Item Number CA138
Item Calcium Chloride, Anhydrous, Granular, FCC
CAS Number 10043-52-4
Molecular Formula CaCl_2
Molecular Weight 110.98
Assay (CaCl_2).....93.0-100.5%

MAXIMUM LIMITS

Acid-Insoluble Matter..... 0.02%
Arsenic (as As)..... 3 mg/kg
Fluoride..... 0.004%
Lead..... 5 mg/kg
Magnesium and Alkali Salts..... 5.0%

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	