



Item Number	CA170
Item	Calcium Sulfate, Dihydrate, NF
CAS Number	10101-41-4
Molecular Formula	$\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$
Molecular Weight	172.17
	Assay (Dried Basis)..... 98.0-101.0%
	Loss on Drying..... 19.0-23.0%
	MAXIMUM LIMITS
	Iron (Fe).....0.01%
	Heavy Metals.....0.001%
	Residual Solvents.....To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843