



Item Number CA194

Item Carboxymethylcellulose Sodium, High

CAS Number 9004-32-4

Molecular Formula

Molecular Weight

Actual Viscosity of 1% Aqueous Solution @ 25°C..... 1,500-3,000 cP

Sodium (Dried Basis)..... 6.5-9.5%

pH of a 1% Solution @ 25°C..... 6.5-8.5

MAXIMUM LIMITS

Loss on Drying..... 10.0%

Heavy Metals..... 0.002%

Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843