



Item Number D1110
Item Dipropylene Glycol Methyl Ether
CAS Number 34590-94-8
Molecular Formula C₇H₁₆O₃
Molecular Weight 148.20

Assay.....Min. 99%
Distillation Range..... 184°-193°C
Specific Gravity @ 25°/25°C..... 0.945-0.957
MAXIMUM LIMITS
Color (APHA)..... 15
Acidity (as CH₃COOH)..... 0.01%
Water..... 0.15%

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	