



Item Number D1146
Item D&C Yellow No. 8
CAS Number 518-47-8
Molecular Formula C₂₀H₁₀O₅Na₂
Molecular Weight 376.28

Assay..... Min. 85.0%

MAXIMUM LIMITS

Sum of Water and Chlorides and Sulfates..... 15%

Matter Insoluble in Alkaline Water.....0.3%

Resorcinol.....0.5%

Phthalic Acid..... 1%

2-(2,4-Dihydroxybenzoyl)benzoic Acid..... 0.5%

Lead (Pb)..... 20 ppm

Arsenic (As).....3 ppm

Mercury (Hg)..... 1 ppm

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843