



Item Number D2440
Item 3,4-Dimethoxybenzaldehyde
CAS Number 120-14-9
Molecular Formula $C_9H_{10}O_3$
Molecular Weight 166.18
Assay.....Approx. 98%
Melting Point.....42°-45°C

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843