



Item Number ET108
Item Alcohol, 190 Proof, USP
CAS Number 64-17-5
Molecular Formula C₂H₅OH
Molecular Weight 46.07

Assay (by Volume)..... 94.9-96.0%
Clarity of Solution..... To pass test
Color of Solution..... To pass test
Specific Gravity @ 15.56°C..... 0.812-0.816
Acidity or Alkalinity..... To pass test

MAXIMUM LIMITS

Ultraviolet Absorption..... To pass test
Organic Impurities..... To pass test
Nonvolatile Residue..... To pass test
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843