



Item Number F1000

Item Ferric Ammonium Citrate, Brown, Powder,

CAS Number 1185-57-5

Molecular Formula

Molecular Weight

Assay (Fe)..... 16.5-18.5%

MAXIMUM LIMITS

Ferric Citrate..... To pass test

Lead..... 2 mg/kg

Mercury..... 1 mg/kg

Oxalate..... To pass test

Sulfate..... 0.3%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843