



Item Number F1015
Item Ferric Chloride, Anhydrous, Sublimed, Powder,
CAS Number 7705-08-0
Molecular Formula FeCl_3
Molecular Weight 162.21

Assay..... Min. 98%

MAXIMUM LIMITS

Alkalies and Earths (as SO_4)..... 0.3%
Copper (Cu).....0.05%
Ferrous Chloride (FeCl_2)..... 1.5%
Insoluble in HCl..... 1.0%
Nitrate (NO_3)..... 0.01%
Phosphorus Compounds (as PO_4).....0.03%
Sulfate (SO_4).....0.01%
Zinc (Zn)..... 0.05%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843