



Item Number IS160
Item Isopropyl Palmitate, NF
CAS Number 142-91-6
Molecular Formula C₁₉H₃₈O₂
Molecular Weight 298.51

Assay..... Min. 90.0%
Specific Gravity @ 25°C.....0.850-0.855
Refractive Index @ 25°C..... 1.435-1.438
Saponification Value..... 183-193

MAXIMUM LIMITS

Acid Value..... 1
Iodine Value..... 1
Residue on Ignition..... 0.1%
Residual Solvents.....To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843