



Item Number L1300

Item L-Lactate Dehydrogenase, Type II, Suspension

CAS Number 9001-60-9

Molecular Formula

Molecular Weight

Activity..... 700-1,200 units/mg protein

MAXIMUM LIMITS

Pyruvate Kinase.....0.01%

Myokinase..... 0.01%

Malic Dehydrogenase.....0.01%

Glutamic-Pyruvic Transaminase..... 0.01%

Glutamic-Oxalacetic Transaminase..... 0.01%

α -Glycerophosphate Dehydrogenase..... 0.01%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843