



Item Number M1023
Item Magnesium Bromide, Hexahydrate, Crystal,
CAS Number 13446-53-2
Molecular Formula $\text{MgBr}_2 \cdot 6\text{H}_2\text{O}$
Molecular Weight 292.25

MAXIMUM LIMITS

Chloride (Cl)..... 0.2%
Iron (Fe).....0.003%
Sulfate (SO_4).....0.05%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843