



Item Number M1114
Item Manganese Sulfate, Monohydrate, Powder,
CAS Number 10034-96-5
Molecular Formula $\text{MnSO}_4 \cdot \text{H}_2\text{O}$
Molecular Weight 169.02

Assay.....98.0-102.0%
Loss on Ignition..... 10.0-13.0%

MAXIMUM LIMITS

Substances not Pptd. by Ammonium Sulfide.....0.5%
Residual Solvents.....To pass test

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	